Arena Worm Tracker

Overview:

We developed a microfluidics-based system for the analysis of *C. elegans* locomotory behavior in response to defined spatial and temporal stimulus patterns. We designed 2 x 2 cm structured arenas with regulated fluid flow that allow *C. elegans* to perform crawling-like locomotion resembling normal behavior on agar surfaces. Stimuli are delivered in the fluid phase and odorevoked behaviors in the arena are automatically captured, segmented and analyzed using the MATLAB scripts available here.

Copyright Notice & Disclaimer:

This source code (MATLAB) is made freely available for research purposes only. The copyright is retained by the authors. Copyright © 2011, Dirk Albrecht. Cori Bargmann. All rights reserved. The use of this software is subject to the following conditions:

- (1) The package and portions thereof may not be sold. Users may modify the software for their own use, but may not redistribute any version of the software other than the original version without the express permission of the author.
- (2) This software package is provided on an "as is" basis. The author in no way warrants either this software or results it may produce. Use of this software is at your own risk.
- (3) The author is under no obligation to provide any services by way of maintenance, updates or corrections for this software.
- (4) Reports or publications resulting from use of this software package must contain an acknowledgment in the form commonly used in academic research.

Any updates to this software will be available at: http://arenawormtracker.dyndns.org

Contact:


dirk.albrecht@gmail.com

Installation Instructions: (assumes familiarity with MATLAB environment)

- 1. Unzip and extract MATLAB library folder ("ArenaWormTracker") to an appropriate location.
- 2. Set MATLAB path (File → Set Path... → Add folder with Subfolders...) to the chosen library location.

MATLAB function and data file notes:

Analysis flowchart:



MATLAB scripts (m-files) accomplish the following basic tasks:

TrackingScript: calls the following 4 functions in sequence:

- 1. **ArenaTracker**: Segment video file to obtain worm tracks. Code adapted from the parallel worm tracker (http://wormsense.stanford.edu/tracker/).
- 2. **SegmentTracks**: Identifies instantaneous behavioral states from worm tracks and morphological data.
- 3. **Ethogram**: Displays ethogram and summarizes behavioral state probability and speed over time. Allows time-shifting to align stimulus pulses.
- 4. **WormDensity**: Summarizes behavioral state and speed data over space and time.

Video tracking, behavior segmentation, and data analysis were performed using MATLAB v7.0.1 with the Image Processing and Statistics toolboxes. Scripts have been tested on Windows XP and 7.

Input files:

- 1. Video file (*.avi)
 - uncompressed or compressed with MATLAB-compatible codec, e.g. Indeo 5
- 2. Experiment information file (*.txt) optional
 - contains information about the experiment. Especially useful for pulse assays to specify the pulse timing and to automatically determine flow properties.

Pulse assay example:

```
Genotype, N2
Stimulus, 1,4,IAA,0.92
, 7,7.5,IAA,0.92
, 8,8.5,IAA,0.92
, 11,12,IAA,0.92
, 13,14,IAA,0.92

NumCycles, 4
CycleLen, 15
FlowVel, 0.88
Delay, 3.53
```

- specifies four 15 min cycles with 0.92µM IAA stimulus from 1-4 min, 7-7.5 min, etc.
- flow velocity = 0.88 mm/s
- time delay between valve switch and stimulus switch upstream of device = 3.53 s

Stripe assay example:

```
Genotype, N2
Stimulus, 0,80,IAA,0,0.92,0,1.84,0
, 80,120,IAA,0
```

- specifies five stripes from 0-80 min, containing from top: 0 μ M, 0.92 μ M, 0 μ M, 1.84 μ M, 0 μ M IAA and 0 μ M IAA throughout from 80-120 min

Output data files:

Gray indicates variables duplicated in previous files. Numbers preceding a period define column numbers. Numbers preceding a colon indicate the meaning of a data or index value.

1. Tracks file (*.mat): Worm position and morphology only

```
AllData
                    Data matrix containing tracking information about every thresholded object, regardless of size
                               1. Frame #
                               2. Object size (pix^2)
                               3-4. Centroid X-Y
                               5. Object eccentricity
DyeData
                    Data matrix for analysis of dye patterns
ExpData
                    Experiment data
File
                    Filenames from ArenaTracker
Tracks.
                    Track structure
          Path
                               x-y position matrix of centroid (pix)
                               vector of frame numbers
          Frames
          Size
                               vector of object size (pix^2)
          Eccentricity
                               vector of eccentricity
                               vector of major axis (pix)
          MajorAxes
                               vector of minor axis (pix)
          MinorAxes
          Orientation
                               vector of orientation (deg)
                               x-y position matrix of bounding box
          Box
          Frame.Image
                               binary image of thresholded object
background
                    background image
```

2. Segmented tracks file (*_seg.mat): Tracks file with worm behavior

Data matrix containing tracking information about every thresholded object, regardless of size AllData 1. Frame # 2. Object size (pix^2) 3-4. Centroid X-Y 5. Object eccentricity Data matrix for analysis of dye patterns DyeData ExpData Experiment data File Filenames from ArenaTracker SegSettings Segmentation settings Tracks. Track structure Path x-y position matrix of centroid (pix) Frames vector of frame numbers Size vector of object size (pix^2) Eccentricity vector of eccentricity vector of major axis (pix) MajorAxes MinorAxes vector of minor axis (pix) Orientation vector of orientation (deg) x-y position matrix of bounding box Box Frame.Image binary image of thresholded object Code segment code (1-10) Distance vector of distance traveled since previous frame (pix) Speed vector of speed (mm/s) vector of path angle (deg) PathAngle PathAngVel vector of path angular velocity (deg/fr) vector of reversing frames Reverse OriginalDistance vector of distance traveled, before pause analysis structure of segment analysis Segment vector of pausing frames Stall NoStall vector of non-pausing frames Beh behavior code 1: Forward straight 2: Forward curve 3: Reverse 4: Pause 5: Pirouette Reverse 6: Pirouette Forward 7: Unknown 8: Out of bounds OriginalTrack original track linked to current track, aka animal # background background image trackstats Tracking statistics 1. Length of track (frames) 2. Animal #

3. Behavior tracks file (*_beh.mat): Condensed version of Segmented tracks file

Tra	acks Track s	structure
	Frames	vector of frame numbers
	Speed	vector of speed (mm/s)
	PathAngle	vector of path angle (deg)
	Beh	behavior code
		1: Forward straight
		2: Forward curve
		3: Reverse
		4: Pause
		5: Pirouette Reverse
		6: Pirouette Forward
		7: Unknown
		8: Out of bounds
	OriginalTrack	original track linked to current track, aka animal #
	X, Ÿ	x- and y-position of centroid

4. Ethogram file (*_ethogram.mat): Worm behavior, time-shifted for pulse assay

```
ExpData
                    Experiment data
                    Data structure containing matrices (row = animal #; col = timepoint) of:
Data.__
         xmat
                    x-position (pix)
                    y-position (pix)
         ymat
                    direction (deg)
         dirmat
         behmat
                   behavior code
                              1: Forward straight
                              2: Forward curve
                              3: Reverse
                              4: Pause
                              5: Pirouette Reverse
                              6: Pirouette Forward
                              7: Unknown
                              8: Out of bounds
         spdmat
                    Speed (mm/s), negative values = reverse.
                    matrix of instantaneous behvaior probability (behavior codes 1-6 only)
                    number of animal-frames ber time point (useful for weighted averages)
         behnum
                    Average instantaneous speed for:
Data.speed.
                              all behaviors
                              forward and pause behaviors
         .fwdpause
         .forward
                              forward runs only
```

5. WormDensity file (*_density.mat): Worm behavior, averaged over time and space

All.	# animal-frames per bin of All animals, including out of bounds and collisions.
	For collisions, number of frames equals estimated number of animals in the collision.
Track	# animal-frames per bin of all valid tracked animals, excluding out-of-bounds,
	collisions, and unknown behaviors.
Beh(1:6)	# animal-frames per bin of animals in a specific behavior:
	1: Forward Run
	2: Pause
	3: Reversal (non-omega)
	4: Pirouette Reverse
	5: Pirouette Forward
	6: Pirouette Rev or Fwd
Speed.FP	
.wt	Weight per bin, i.e. number of animal-frames contributing to each speed average per bin
XYTimeD	4-D matrix, Xbin, Ybin, Timebin, Direction (1=up, 2=horiz, 3=down)
XYTime	3-D matrix, sum of XYTimeD over all D
XY	2-D
TimeX	2-D
TimeY	2-D
XYD	3-D
XTimeD	3-D
YTimeD	3-D
Time	1-D vector of time bins (mins)
X Y	1-D vector of x-position (mm) 1-D vector of y-position (mm)
. I	1-b vector or y-position (min)
ExpData	Experiment data
File	Filenames from ArenaTracker
FileInfo	Additional filenames
SegData	Data matrix of all segments per experiment
	1. Track #
	2. Segment # (of Track)
	3. Behavior Code
	4. Behavior Run (Run = 1 vs. Turn = 0)
0 5	5. Segment Length (frames)
SegPos	Position matrix of all segments per experiment
	1-3. SegmentStart (frame, X, Y of starting point)

	4-6. SegmentEnd (frame, X, Y of ending point)
	7-8. SegmentCentroid (X, Y of centroid)
	9-12. SegmentMinMax (minX, minY, maxX, maxY)
SegDir	Direction matrix of all segments per experiment
	1. angUD (ending direction of segment: 1=up (-y); 2=horiz; 3=dn (+y)
	2-3. SegmentDir (mean, std of all angles in segment: +x axis; ccw)
	4. SegmentBoundHead (initial, final angle)
	5-17. Histogram of angles (# frames in 0,30,60,90, 330)
TrackData	Data matrix of every animal and time point during the experiment
	1. Animal#
	2. Time (min)
	3-4. XY centroid position
	5. Animal Behavior
	1: out of bounds
	2: forward run
	3: pause
	4: reversal
	5: omega reverse
	6: omega forward
	6. Speed (forward & pause only, otherwise NaN)
	7. Animal Direction Category
	-1: -y ("up")
	0: horizontal
	1: +y
	2: no direction
	8. Animal Direction angle (deg)
TurnData	Data matrix of every turn during the experiment
	1. time(fr)
	2. posX
	3. posY
	4. turnID
	1: forward run (2= with reverse frames)
	3: forward shift
	4: 60deg curve
	5: unknown short turning segment
	6: long curve (non-pirouette)
	7: pirouette
	8: unknown long turning segment
	5. priorDirUD
	6. afterDirUD
segidx	segment # of forward run prior to turn for each row in TurnData