

Actual observations:

Without MB

With MB

Both

**I. HARVEST** (Gandre, S., & Bliet, A. M., 2016)

Materials

~~5 worm plates per harvest~~ [see next bullet]

- 1 (large) plate per treatment group

- 0.1M NaCl

- M9 buffer

~~Mitochondrial IB~~

~~210 mM mannitol~~

~~70 mM sucrose~~

~~0.1 mM EDTA, pH 8.0~~

~~5 mM Tris HCl, pH 7.4~~

~~1 mM PMSF~~

- Centrifuge

1. Harvest worms from 5 plates (either MB treated or untreated) by washing off the plate with 10 mL M9 buffer per plate

a. x3 10mL/plate + 3mL to catch stragglers

2. Pellet the worms (all 5 plates together) by centrifugation for 5 min @ ~~1000g~~ 2500g

a. Loose pellet formed; sent in for +10 min @ 2500g

b. Pellet still loose (especially blue); + 10 min @ 5000g

c. + 10 min @ 5000g

3. Resuspend the pellet in 10 volumes ice-cold 0.1M NaCl

a. Add an equal volume of ice-cold 60% sucrose solution; mix by inversion

b. Split into x2 1 mL tubes/group

4. Centrifuge for 2 min @ 500g

5. Transfer the uppermost layer (floating worms) to a new tube using a Pasteur pipet

a. Pool top 250µL from each tube (each group) → 1 tube per group

~~6. Dilute the worm suspension with 10 volumes ice-cold 0.1M NaCl~~

a. Centrifuge for 2.5 min @ ~~500g~~ 5000g

b. +10 mins @ 10000g

~~7. Discard the supernatant and resuspend the pellet in 50 volumes ice-cold 0.1M NaCl~~

**II. ISOLATE** (Gandre, S., & Bliet, A. M., 2016)

Materials (volumes were adjusted to reduced sample quantity)

- 5 plates worth of worms, pelleted

- Mitochondrial Isolation Buffer (IB) (for 25 mL)

o 210 mM mannitol -- 0.956 g

o 70 mM sucrose -- 0.599 g

- 0.1 mM EDTA, pH 8.0 -- 0.00093 g (0.93 mg)
  - 5 mM Tris-HCl, pH 7.4 -- 0.0151 g (15.1 mg)
  - 1 mM PMSF -- 0.00435 g (4.35 mg) -- fume hood?
  - Potter-Elvehjem homogenizer/Dounce/honestly, I will use a rock at this point
  - Centrifuge
1. Pellet worms by centrifuging for 5 min @ 2000g
  2. Resuspend 5 g worms (this was originally done for twice the amount of sample, so it's probably going to be *way* less than 5 g) in 5 mL ice-cold IB
  3. Homogenize with 15 strokes of a chilled Potter-Elvehjem homogenizer
  4. Collect the homogenate in a 25 mL Falcon tube and increase the volume to 13 mL with IB
    - a. Centrifuge for 10 min @ 750g
  5. Transfer the supernatant to a fresh tube and resuspend the pellet in 5 mL IB
  6. Homogenize the remaining pellet with 15 more strokes of the homogenizer
    - a. Repeat Steps 4 & 5
  7. Combine the supernatants from Steps 5 & 6a; ~~save an aliquot as total worm lysate~~
  8. Centrifuge the supernatants for 10 min @ 12,000g; ~~save an aliquot of the resulting supernatant as postmitochondrial supernatant~~
  9. Resuspend the mitochondrial pellet in 6 mL IB
    - a. Centrifuge for 10 min @ 750g
  10. Transfer the supernatant to a new tube without disrupting the pellet
    - a. Centrifuge for 10 min @ 12,000g and discard the supernatant
  11. Combine and resuspend the mitochondrial pellets from Steps 9a & 10a in IB

### III. ASSAY (Sigma-Aldrich, 2019)

#### Materials

- Isolated *C. elegans* mitochondria (MB treated or untreated)
  - Equine heart cytochrome c
  - Reaction buffer, on ice
    - 50 mM Tris-HCl, pH 7.4 -- 6.057 g
    - 250 mM sucrose -- 85.575 g
    - 10 mM KCl -- 0.746 g
    - 1 mM EDTA -- 0.372 g
    - 10 mM KH<sub>2</sub>PO<sub>4</sub> -- 1.361 g
    - Final volume 1L w/distilled water
    - Filter-sterilize through a 0.22 μm filter (optional)
  - Spectrophotometer
1. Add ~~50 μL 50 μM~~ 50 μL 4.4 mM oxidized cytochrome c to 950 μL of buffer in a cuvette (blank)
    - a.
  2. Add ~~50 μM~~ 50 μL 4.4 mM oxidized cytochrome c to 950 μL of buffer to a new cuvette
    - a. Add 50 μL mitochondria
  3. Allow the reactions to proceed at 25°C for 5 minutes; monitor absorbance

4. Add dithiothreitol (DTT) to a cuvette containing 950  $\mu\text{L}$  of buffer and 50  $\mu\text{L}$  4.4 mM oxidized cytochrome c to observe the full reduction of cytochrome c.

5. Calculate reduction rate

a. 
$$\frac{\Delta A}{\text{min}} = (A/\text{min})_{\text{sample}} - (A/\text{min})_{\text{blank}}$$

b. 
$$\text{Enzyme Activity (Units/mL)} = \frac{\Delta A/\text{min} * \text{dil} * V_{\text{reaction}}}{V_{\text{enzyme}} * 21.84}$$

where:

21.84 = extinction coefficient of reduced cytochrome c ( $\text{mM}^{-1}\text{cm}^{-1}$ )

$V_{\text{reaction}}$  = total reaction volume

$V_{\text{enzyme}}$  = volume of mitochondrial sample

dil = dilution of mitochondrial sample

6. Compare treatment groups

a. 
$$\text{Difference (\%)} = \frac{\text{Activity}_{\text{MB}} - \text{Activity}_{\text{untreated}}}{\text{Activity}_{\text{untreated}}} * 100$$